

For Research Use Only

# NEGR1 Polyclonal antibody

Catalog Number: 13674-1-AP

2 Publications



## Basic Information

Catalog Number:

13674-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4659

GenBank Accession Number:

BC036771

GeneID (NCBI):

257194

UNIPROT ID:

Q7Z3B1

Full Name:

neuronal growth regulator 1

Calculated MW:

226 aa, 25 kDa

Observed MW:

46 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : mouse cerebellum tissue, rat cerebellum tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yu-Han H Hsu	37207277	iScience	
Pasquale Nigro	37058410	Cell Rep	IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 13674-1-AP (NEGR1 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.