

For Research Use Only

# Glypican 5 Polyclonal antibody

Catalog Number: 13678-1-AP



## Basic Information

Catalog Number:

13678-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4493

GenBank Accession Number:

BC039730

GeneID (NCBI):

2262

UNIPROT ID:

P78333

Full Name:

glypican 5

Calculated MW:

572 aa, 64 kDa

Observed MW:

64 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

## Background Information

GPC5, also named as Glypican-5 and Secreted glypican-5, is cell surface proteoglycan that bears heparan sulfate. It is primarily expressed in the brain. 65kd band is Glycosylation; 45kd band is mature form.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

