

For Research Use Only

PRKACG Polyclonal antibody

Catalog Number: 13697-1-AP



Basic Information

Catalog Number:

13697-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4595

GenBank Accession Number:

BC039888

GeneID (NCBI):

5568

UNIPROT ID:

P22612

Full Name:

protein kinase, cAMP-dependent,
catalytic, gamma

Calculated MW:

351 aa, 40 kDa

Observed MW:

38 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: human testis tissue, mouse testis tissue, rat testis
tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 13697-1-AP (PRKACG antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.