

For Research Use Only

NDST2 Polyclonal antibody

Catalog Number: 13716-1-AP



Basic Information

Catalog Number: 13716-1-AP	GenBank Accession Number: BC035711	Purification Method: Antigen affinity purification
Size: 190 µg/ml	GeneID (NCBI): 8509	
Source: Rabbit	UNIPROT ID: P52849	
Isotype: IgG	Full Name: N-deacetylase/N-sulfotransferase (heparan glucosaminyl) 2	
Immunogen Catalog Number: AG4656	Calculated MW: 883 aa, 101 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data