

For Research Use Only

RAX2 Polyclonal antibody

Catalog Number: 13735-1-AP



Basic Information

Catalog Number:

13735-1-AP

Size:

100 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4589

GenBank Accession Number:

BC032512

GeneID (NCBI):

84839

UNIPROT ID:

Q96IS3

Full Name:

retina and anterior neural fold
homeobox 2

Calculated MW:

184 aa, 20 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

