

For Research Use Only

KCNIP4 Polyclonal antibody

Catalog Number: 13748-1-AP



Basic Information

Catalog Number:

13748-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4688

GenBank Accession Number:

BC032520

GeneID (NCBI):

80333

UNIPROT ID:

Q6PIL6

Full Name:

Kv channel interacting protein 4

Calculated MW:

250 aa, 28 kDa

Observed MW:

29 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : rat brain tissue, mouse brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat brain tissue were subjected to SDS PAGE followed by western blot with 13748-1-AP (KCNIP4 Antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.