

For Research Use Only

RAB15 Polyclonal antibody

Catalog Number: 13800-1-AP

2 Publications



Basic Information

Catalog Number:

13800-1-AP

Size:

160 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4775

GenBank Accession Number:

BC040679

GeneID (NCBI):

376267

UNIPROT ID:

P59190

Full Name:

RAB15, member RAS oncogene family

Calculated MW:

208 aa, 24 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

mouse, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Qingli Liu	36309191	Acta Biomater	IHC
Xiao-Hu Wang	35286985	Arch Gerontol Geriatr	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

