

For Research Use Only

CCL21 Polyclonal antibody, PBS Only

Catalog Number:13814-1-PBS



Basic Information

Catalog Number: 13814-1-PBS	GenBank Accession Number: BC027918	Purification Method: Antigen affinity purification
Size: 1 mg/ml	GeneID (NCBI): 6366	
Source: Rabbit	UNIPROT ID: O00585	
Isotype: IgG	Full Name: chemokine (C-C motif) ligand 21	
Immunogen Catalog Number: AG4897	Calculated MW: 134 aa, 15 kDa	

Applications

Tested Applications:
Indirect ELISA
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data