

For Research Use Only

CYTSB Polyclonal antibody

Catalog Number: 13851-1-AP



Basic Information

Catalog Number:

13851-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4822

GenBank Accession Number:

BC050058

GeneID (NCBI):

92521

UNIPROT ID:

Q5M775

Full Name:

cytospin B

Calculated MW:

119 kDa

Observed MW:

119 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse kidney tissue, HeLa cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 13851-1-AP (CYTSB antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.