

For Research Use Only

Calpain 7 Polyclonal antibody

Catalog Number: 13870-1-AP



Basic Information

Catalog Number:

13870-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4820

GenBank Accession Number:

BC056202

GeneID (NCBI):

23473

UNIPROT ID:

Q9Y6W3

Full Name:

calpain 7

Calculated MW:

93 kDa

Observed MW:

92 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human heart tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human heart tissue were subjected to SDS PAGE followed by western blot with 13870-1-AP (Calpain 7 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.