

For Research Use Only

ZHX1 Polyclonal antibody, PBS Only

Catalog Number: 13903-1-PBS

Featured Product



Basic Information

Catalog Number:

13903-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4874

GenBank Accession Number:

BC040481

GeneID (NCBI):

11244

UNIPROT ID:

Q9UKY1

Full Name:

zinc fingers and homeoboxes 1

Calculated MW:

873 aa, 98 kDa

Observed MW:

100-120 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Binding of the nuclear factor-Y complex (NF-Y) to the inverted CCAAT-box interferes with transcription activation through nucleosome reorganization. The three homologous proteins forming the zinc-fingers and homeoboxes (ZHX) family interact with the activation domain of NF-Ya to repress transcription [PMID:19348505,20509910]. Human zinc-fingers and homeoboxes 1 (ZHX1) was cloned as a protein that interacts with the activation domain (AD) of the A subunit of nuclear factor-Y (NF-YA) [PMID:12237128].

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 13903-1-AP (ZHX1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 13903-1-PBS in a different storage buffer formulation.



IP result of anti-ZHX1 (IP:13903-1-AP, 4ug; Detection:13903-1-AP 1:500) with HEK-293 cells lysate 3200ug. This data was developed using the same antibody clone with 13903-1-PBS in a different storage buffer formulation.