

For Research Use Only

SLAIN2 Polyclonal antibody

Catalog Number: 13908-1-AP **1 Publications**



Basic Information

Catalog Number:

13908-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4888

GenBank Accession Number:

BC040993

GeneID (NCBI):

57606

UNIPROT ID:

Q9P270

Full Name:

SLAIN motif family, member 2

Calculated MW:

581 aa, 63 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : mouse brain tissue, Transfected HEK-293 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Elif Nur Firat-Karalar	32256142	Turk J Biol	WB, IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 13908-1-AP (SLAIN2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.