

For Research Use Only

# JIP2 Polyclonal antibody, PBS Only

Catalog Number: 13912-1-PBS



## Basic Information

Catalog Number:

13912-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4902

GenBank Accession Number:

BC047527

GeneID (NCBI):

23542

UNIPROT ID:

Q13387

Full Name:

mitogen-activated protein kinase 8  
interacting protein 2

Calculated MW:

824 aa, 88 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data