

For Research Use Only

ZnT7 Polyclonal antibody

Catalog Number: 13966-1-AP **9 Publications**



Basic Information

Catalog Number: 13966-1-AP	GenBank Accession Number: BC064692	Purification Method: Antigen affinity purification
Size: 300 µg/ml	GeneID (NCBI): 148867	Recommended Dilutions: WB 1:500-1:2000 IHC 1:50-1:500 IF 1:20-1:200
Source: Rabbit	UNIPROT ID: Q8NEW0	
Isotype: IgG	Full Name: solute carrier family 30 (zinc transporter), member 7	
Immunogen Catalog Number: AG4957	Calculated MW: 42 kDa	
	Observed MW: 42 kDa	

Applications

Tested Applications: IF/ICC, IHC, WB, ELISA	Positive Controls: WB : C2C12 cells, IHC : human lung cancer tissue, IF : A431 cells,
Cited Applications: WB, IF, IHC	
Species Specificity: human, mouse, rat	
Cited Species: chicken, human, mouse, rat	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Masayuki Kaneko	25284286	J Neurosci Res	WB, IF
Qian Sun	25767260	Am J Physiol Gastrointest Liver Physiol	WB
Surapun Tepasomdech	28469980	FEBS Open Bio	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 13966-1-AP (ZnT7 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 13966-1-AP (ZnT7 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 13966-1-AP (ZnT7 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of A431 cells, using SLC30A7 antibody 13966-1-AP at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).