

For Research Use Only

Calicin Polyclonal antibody

Catalog Number: 13972-1-AP

1 Publications



Basic Information

Catalog Number:

13972-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5005

GenBank Accession Number:

BC019251

GeneID (NCBI):

881

UNIPROT ID:

Q13939

Full Name:

calicin

Calculated MW:

67 kDa

Observed MW:

60-65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : human testis tissue, rat testis tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jiaxin He	36527329	Clin Genet	WB, IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 13972-1-AP (Calicin antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.