

For Research Use Only

FZD6 Polyclonal antibody

Catalog Number: 13982-1-AP



Basic Information

Catalog Number: 13982-1-AP	GenBank Accession Number: BC060836	Purification Method: Antigen affinity purification
Size: 190 µg/ml	GeneID (NCBI): 8323	
Source: Rabbit	UNIPROT ID: O60353	
Isotype: IgG	Full Name: frizzled homolog 6 (Drosophila)	
Immunogen Catalog Number: AG5027	Calculated MW: 79 kDa	

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human breast cancer using 13982-1-AP (FZD6 antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human breast cancer using 13982-1-AP (FZD6 antibody) at dilution of 1:50 (under 40x lens).