

For Research Use Only

# B3GNT4 Polyclonal antibody

Catalog Number: 13984-1-AP



## Basic Information

Catalog Number:

13984-1-AP

Size:

170 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5050

GenBank Accession Number:

BC031399

GeneID (NCBI):

79369

UNIPROT ID:

Q9C0J1

Full Name:

UDP-GlcNAc:betaGal beta-1,3-N-acetylglucosaminyltransferase 4

Calculated MW:

378 aa, 42 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

