

For Research Use Only

GOLPH3 Polyclonal antibody

Catalog Number: 14053-1-AP



Basic Information

Catalog Number: 14053-1-AP	GenBank Accession Number: BC033725	Purification Method: Antigen affinity purification
Size: 150 µg/ml	GeneID (NCBI): 64083	
Source: Rabbit	UNIPROT ID: Q9H4A6	
Isotype: IgG	Full Name: golgi phosphoprotein 3 (coat-protein)	
Immunogen Catalog Number: AG5278	Calculated MW: 298 aa, 34 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data