

For Research Use Only

DLK Polyclonal antibody

Catalog Number: 14054-1-AP



Basic Information

| | | |
|--|---|--|
| Catalog Number: 14054-1-AP | GenBank Accession Number: BC050050 | Purification Method: Antigen affinity purification |
| Size: 330 µg/ml | GeneID (NCBI): 7786 | |
| Source: Rabbit | UNIPROT ID: Q12852 | |
| Isotype: IgG | Full Name: mitogen-activated protein kinase kinase kinase 12 | |
| Immunogen Catalog Number: AG4780 | Calculated MW: 567 aa, 64 kDa | |

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data