

For Research Use Only

AMIGO2 Polyclonal antibody

Catalog Number: 14076-1-AP



Basic Information

Catalog Number:

14076-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5200

GenBank Accession Number:

BC047595

GeneID (NCBI):

347902

UNIPROT ID:

Q86S12

Full Name:

adhesion molecule with Ig-like domain 2

Calculated MW:

58 kDa

Observed MW:

58 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:300-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Y79 cells, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Y79 cells were subjected to SDS PAGE followed by western blot with 14076-1-AP (AMIGO2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.