

For Research Use Only

OVOL1 Polyclonal antibody

Catalog Number: 14082-1-AP

Featured Product

7 Publications



Basic Information

Catalog Number:

14082-1-AP

Concentration:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5210

GenBank Accession Number:

BC059408

GeneID (NCBI):

5017

UNIPROT ID:

O14753

Full Name:

ovo-like 1(Drosophila)

Calculated MW:

30 kDa

Observed MW:

30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB : mouse brain tissue, human brain tissue, rat brain tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Hernan Roca	24124593	PLoS One	WB
Lu Chen	29725020	Cell Death Dis	WB
Fernando Eguiarte-Solomon	33713497	Mol Carcinog	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14082-1-AP (OVOL1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 14082-1-AP (OVOL1 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.