

For Research Use Only

IRGC Polyclonal antibody

Catalog Number: 14090-1-AP

1 Publications



Basic Information

Catalog Number:

14090-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5279

GenBank Accession Number:

BC066939

GeneID (NCBI):

56269

UNIPROT ID:

Q6NXR0

Full Name:

immunity-related GTPase family,
cinema

Calculated MW:

50 kDa

Observed MW:

50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : mouse testis tissue, rat testis tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yuki Kaneda	35618043	Dev Biol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14090-1-AP (IRGC antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.