

For Research Use Only

ASB17 Polyclonal antibody

Catalog Number: 14146-1-AP



Basic Information

Catalog Number:

14146-1-AP

Size:

220 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5315

GenBank Accession Number:

BC036573

GeneID (NCBI):

127247

UNIPROT ID:

Q8WVJ9

Full Name:

ankyrin repeat and SOCS box-containing 17

Calculated MW:

295 aa, 34 kDa

Observed MW:

34 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse testis tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue lysates were subjected to SDS PAGE followed by western blot with 14146-1-AP (ASB17 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.