

For Research Use Only

PTGR2 Polyclonal antibody

Catalog Number: 14164-1-AP



Basic Information

Catalog Number:

14164-1-AP

Size:

750 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5352

GenBank Accession Number:

BC059364

GeneID (NCBI):

145482

UNIPROT ID:

Q8N8N7

Full Name:

prostaglandin reductase 2

Calculated MW:

351 aa, 39 kDa

Observed MW:

38.5 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse kidney tissue, BxPC-3 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 14164-1-AP (PTGR2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.