

For Research Use Only

POTEB Polyclonal antibody

Catalog Number: 14169-1-AP

1 Publications



Basic Information

Catalog Number:

14169-1-AP

Size:

110 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5368

GenBank Accession Number:

BC066981

GeneID (NCBI):

339010

UNIPROT ID:

Q6S5H4

Full Name:

POTE ankyrin domain family,
member B

Calculated MW:

66 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Cited Applications:

IHC

Species Specificity:

human

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Redfield Samantha M SM	23841027	Am J Cancer Res	IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

