

For Research Use Only

SMC5 Polyclonal antibody, PBS Only

Catalog Number: 14178-1-PBS

Featured Product



Basic Information

Catalog Number:

14178-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5385

GenBank Accession Number:

BC038225

GeneID (NCBI):

23137

UNIPROT ID:

Q8IY18

Full Name:

structural maintenance of
chromosomes 5

Calculated MW:

1101 aa, 129 kDa

Observed MW:

120-130 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 14178-1-AP (SMC5 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14178-1-PBS in a different storage buffer formulation.