

For Research Use Only

RBMV1A1 Polyclonal antibody

Catalog Number: 14239-1-AP



Basic Information

Catalog Number:

14239-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5493

GenBank Accession Number:

BC047768

GeneID (NCBI):

5940

UNIPROT ID:

P0DJD3

Full Name:

RNA binding motif protein, Y-linked, family 1, member A1

Calculated MW:

56 kDa

Observed MW:

50-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : COLO 320 cells, NCCIT cell

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



COLO 320 cells were subjected to SDS PAGE followed by western blot with 14239-1-AP (RBMV1A1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.