

For Research Use Only

APLF Polyclonal antibody

Catalog Number: 14252-1-AP

1 Publications



Basic Information

Catalog Number:

14252-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5525

GenBank Accession Number:

BC041144

GeneID (NCBI):

200558

UNIPROT ID:

Q8IW19

Full Name:

apratxin and PNKP like factor

Calculated MW:

57 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : PC-3 cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Matthew C J Yip	31011209	Nat Struct Mol Biol	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 14252-1-AP (APLF antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.