

For Research Use Only

APLF Polyclonal antibody

Catalog Number: 14252-1-AP **1 Publications**



Basic Information

Catalog Number: 14252-1-AP	GenBank Accession Number: BC041144	Purification Method: Antigen affinity purification
Size: 500 µg/ml	GeneID (NCBI): 200558	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q8IW19	
Isotype: IgG	Full Name: apratxin and PNKP like factor	
Immunogen Catalog Number: AG5525	Calculated MW: 57 kDa	
	Observed MW: 60 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : PC-3 cells,
Species Specificity: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Matthew C J Yip	31011209	Nat Struct Mol Biol	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 14252-1-AP (APLF antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.