

For Research Use Only

VPS36 Polyclonal antibody, PBS Only

Catalog Number:14262-1-PBS



Basic Information

Catalog Number:

14262-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5549

GenBank Accession Number:

BC050439

GeneID (NCBI):

51028

UNIPROT ID:

Q86VN1

Full Name:

vacuolar protein sorting 36 homolog
(S. cerevisiae)

Calculated MW:

44 kDa

Observed MW:

43 kDa, 65 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HL-60 cells were subjected to SDS PAGE followed by western blot with 14262-1-AP (VPS36 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14262-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 14262-1-AP (VPS36 antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 14262-1-PBS in a different storage buffer formulation.