

For Research Use Only

FBXO18 Polyclonal antibody

Catalog Number: 14275-1-AP



Basic Information

Catalog Number: 14275-1-AP	GenBank Accession Number: BC032674	Purification Method: Antigen affinity purification
Size: 300 µg/ml	GeneID (NCBI): 84893	Recommended Dilutions: WB 1:500-1:1000 IF 1:50-1:500
Source: Rabbit	UNIPROT ID: Q8NFZO	
Isotype: IgG	Full Name: F-box protein, helicase, 18	
Immunogen Catalog Number: AG5622	Calculated MW: 118 kDa, 123 kDa	
	Observed MW: 118-123 kDa	

Applications

Tested Applications: IF/ICC, WB, ELISA	Positive Controls: WB : PC-3 cells, HeLa cells IF : PC-3 cells,
Species Specificity: human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 14275-1-AP (FBXO18 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (10% Formaldehyde) fixed PC-3 cells using 14275-1-AP (FBXO18 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).