

For Research Use Only

FBXO18 Polyclonal antibody, PBS Only

Catalog Number: 14275-1-PBS



Basic Information

Catalog Number:

14275-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5622

GenBank Accession Number:

BC032674

GeneID (NCBI):

84893

UNIPROT ID:

Q8NFZ0

Full Name:

F-box protein, helicase, 18

Calculated MW:

118 kDa, 123 kDa

Observed MW:

118-123 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 14275-1-AP (FBXO18 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14275-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (10% Formaldehyde) fixed PC-3 cells using 14275-1-AP (FBXO18 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 14275-1-PBS in a different storage buffer formulation.