

For Research Use Only

# SUPT3H Polyclonal antibody

Catalog Number: 14279-1-AP



## Basic Information

Catalog Number:

14279-1-AP

Size:

350 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5434

GenBank Accession Number:

BC050384

GeneID (NCBI):

8464

UNIPROT ID:

O75486

Full Name:

suppressor of Ty 3 homolog (S. cerevisiae)

Calculated MW:

44 kDa

Observed MW:

34-37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : A431 cells, HL-60 cells, MCF-7 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14279-1-AP (SUPT3H antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.