

For Research Use Only

BMPR1B Polyclonal antibody

Catalog Number: 14312-1-AP



Basic Information

Catalog Number:

14312-1-AP

Size:

160 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5545

GenBank Accession Number:

BC047773

GeneID (NCBI):

658

UNIPROT ID:

O00238

Full Name:

bone morphogenetic protein receptor,
type IB

Calculated MW:

57 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

