

For Research Use Only

DVL1 Polyclonal antibody

Catalog Number: 14314-1-AP

1 Publications



Basic Information

Catalog Number:

14314-1-AP

Size:

70 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5550

GenBank Accession Number:

BC050454

GeneID (NCBI):

1855

UNIPROT ID:

O14640

Full Name:

dishevelled, dsh homolog 1
(Drosophila)

Calculated MW:

75 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

rat, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yiyang Du	37739105	J Ethnopharmacol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

