

For Research Use Only

# CLCA2 Polyclonal antibody

Catalog Number: 14323-1-AP

2 Publications



## Basic Information

Catalog Number:

14323-1-AP

Size:

140 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5582

GenBank Accession Number:

BC041096

GeneID (NCBI):

9635

UNIPROT ID:

Q9UQC9

Full Name:

chloride channel accessory 2

Calculated MW:

104 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human, mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Linda D Thomas	32873706	J Biol Chem	WB
Linda D Thomas	37225015	Mol Metab	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

