

For Research Use Only

# BPI Polyclonal antibody

Catalog Number: 14368-1-AP



## Basic Information

Catalog Number:

14368-1-AP

Size:

160 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5695

GenBank Accession Number:

BC040955

GeneID (NCBI):

671

Full Name:

bactericidal/permeability-increasing protein

Calculated MW:

54 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF 1:10-1:100

## Applications

Tested Applications:

IF, ELISA

Species Specificity:

human

Positive Controls:

IF : HepG2 cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of HepG2 cells, using BPI antibody 14368-1-AP at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).