

For Research Use Only

BPI Polyclonal antibody, PBS Only

Catalog Number:14368-1-PBS



Basic Information

Catalog Number:

14368-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5695

GenBank Accession Number:

BC040955

GeneID (NCBI):

671

UNIPROT ID:

P17213

Full Name:

bactericidal/permeability-increasing protein

Calculated MW:

54 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 14368-1-AP (BPI antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 14368-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HL-60 cells using BPI antibody (14368-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 14368-1-PBS in a different storage buffer formulation.