

For Research Use Only

NLK Polyclonal antibody

Catalog Number: 14374-1-AP

Featured Product

3 Publications



Basic Information

Catalog Number:

14374-1-AP

Size:

210 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5750

GenBank Accession Number:

BC064663

GeneID (NCBI):

51701

UNIPROT ID:

Q9UBE8

Full Name:

nemo-like kinase

Calculated MW:

58 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Dongling Zou	25833695	Tumour Biol	WB
Guozhen Wang	38910259	J Nanobiotechnology	WB, IHC
Cemre Ural	37805666	Amino Acids	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

