

For Research Use Only

# BIK Polyclonal antibody

Catalog Number: 14389-1-AP



## Basic Information

Catalog Number:

14389-1-AP

Size:

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5859

GenBank Accession Number:

BC001599

GeneID (NCBI):

638

UNIPROT ID:

Q13323

Full Name:

BCL2-interacting killer (apoptosis-inducing)

Calculated MW:

18 kDa

Purification Method:

Antigen Affinity purified

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

BIK, also named as NBK, BIP1 and BP4, is a pro-apoptotic BCL-2 family member and is the founding member of a subfamily of pro-apoptotic proteins known as BH3-alone proteins. It is a transcript that is upregulated during induction of apoptosis.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

