

For Research Use Only

SETDB2 Polyclonal antibody

Catalog Number: 14428-1-AP

2 Publications



Basic Information

Catalog Number:

14428-1-AP

Size:

233 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5795

GenBank Accession Number:

BC047434

GeneID (NCBI):

83852

UNIPROT ID:

Q96T68

Full Name:

SET domain, bifurcated 2

Calculated MW:

82 kDa

Observed MW:

101 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : Jurkat cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yanqing Liu	39631394	Mol Cell	WB
Jingyao Zhang	37197420	Front Oncol	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 14428-1-AP (SETDB2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.