

For Research Use Only

beta-actin-Specific Polyclonal antibody



Catalog Number: 14439-1-AP **5 Publications**

Basic Information

Catalog Number: 14439-1-AP	GenBank Accession Number: NM_001101	Purification Method: Antigen affinity purification
Size: 70 µg/ml	GeneID (NCBI): 60	
Source: Rabbit	UNIPROT ID: P60709	
Isotype: IgG	Full Name: actin, beta	
	Calculated MW: 42 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse, rat, zebrafish

Cited Species:
human

Background Information

This antibody is specific to B-actin, but not other actins.

Notable Publications

Author	Pubmed ID	Journal	Application
P Luo	29740907	Neurogastroenterol Motil	WB
Eman S Abou Zaid	36600460	Int J Immunopathol Pharmacol	WB
Meng Dai	36529788	Eur J Med Res	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data