

For Research Use Only

# beta-actin-Specific Polyclonal antibody



Catalog Number: 14439-1-AP

6 Publications

## Basic Information

Catalog Number:

14439-1-AP

Size:

70 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001101

GeneID (NCBI):

60

UNIPROT ID:

P60709

Full Name:

actin, beta

Calculated MW:

42 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat, zebrafish

Cited Species:

human

## Background Information

This antibody is specific to B-actin, but not other actins.

## Notable Publications

Author	Pubmed ID	Journal	Application
P Luo	29740907	Neurogastroenterol Motil	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

