

For Research Use Only

RELT Polyclonal antibody

Catalog Number: 14532-1-AP



Basic Information

Catalog Number:

14532-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6032

GenBank Accession Number:

BC051810

GeneID (NCBI):

84957

UNIPROT ID:

Q969Z4

Full Name:

RELT tumor necrosis factor receptor

Calculated MW:

46 kDa

Observed MW:

46-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2400

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : PC-3 cells, HeLa cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 14532-1-AP (RELT antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.