

For Research Use Only

RELT Polyclonal antibody

Catalog Number: 14532-1-AP



Basic Information

Catalog Number:

14532-1-AP

Size:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6032

GenBank Accession Number:

BC051810

GeneID (NCBI):

84957

UNIPROT ID:

Q969Z4

Full Name:

RELT tumor necrosis factor receptor

Calculated MW:

46 kDa

Observed MW:

46 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Daudi cells, HeLa cells, HL-60 cells, Jurkat cells, Raji cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14532-1-AP (RELT antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.