

For Research Use Only

CASZ1 Polyclonal antibody, PBS Only

Catalog Number:14534-1-PBS



Basic Information

Catalog Number: 14534-1-PBS	GenBank Accession Number: BC051883	Purification Method: Antigen affinity purification
Source: Rabbit	GeneID (NCBI): 54897	
Isotype: IgG	UNIPROT ID: Q86V15	
Immunogen Catalog Number: AG6035	Full Name: castor zinc finger 1	
	Calculated MW: 190 kDa	
	Observed MW: 150 kDa, 250 kDa	

Applications

Tested Applications:
WB, IP, Indirect ELISA

Species Specificity:
human, mouse

Background Information

Zinc finger protein castor homolog 1 (CASZ1), also known as zinc finger protein 693 (ZNF693), belongs to the C2H2-type zinc finger protein family. CASZ1 is a transcription factor that is evolutionarily conserved and participates in multiple embryonic development and physiological processes. CASZ1 has been found to regulate aldosterone synthesis and repress the activation of aldosterone target genes, suggesting that it may participate in blood pressure regulation. CASZ1 is ubiquitously expressed in various embryonic tissues. Subcellularly, CASZ1 is usually localized in the cell nucleus, consistent with its role in transcription regulation (PMID: 37509718). CASZ1 has two isoforms namely CASZ1a (i.e., hcasz11, 190 kDa) and CASZ1b (i.e., hcasz5, 125 kDa). The observed molecular weights were 150 kDa and 250 kDa, respectively, consistent with descriptions in the literature (PMID: 32060262).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14534-1-AP (CASZ1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14534-1-PBS in a different storage buffer formulation.



IP result of anti-CASZ1 (IP:14534-1-AP, 4ug; Detection:14534-1-AP 1:1000) with C2C12 cells lysate 1200 ug. This data was developed using the same antibody clone with 14534-1-PBS in a different storage buffer formulation.



IP result of anti-CASZ1 (IP:14534-1-AP, 4ug; Detection:14534-1-AP 1:1000) with HUVEC cells lysate 1080 ug. This data was developed using the same antibody clone with 14534-1-PBS in a different storage buffer formulation.