

For Research Use Only

NEUROD4 Polyclonal antibody

Catalog Number: 14610-1-AP



Basic Information

Catalog Number:

14610-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6183

GenBank Accession Number:

BC040961

GeneID (NCBI):

58158

UNIPROT ID:

Q9HD90

Full Name:

neurogenic differentiation 4

Calculated MW:

37 kDa

Observed MW:

50-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Y79 cells, Neuro-2a cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Y79 cells were subjected to SDS PAGE followed by western blot with 14610-1-AP (NEUROD4 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.