

For Research Use Only

# TMEM54 Polyclonal antibody

Catalog Number: 14618-1-AP



## Basic Information

Catalog Number:

14618-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6234

GenBank Accession Number:

BC001418

GeneID (NCBI):

113452

UNIPROT ID:

Q969K7

Full Name:

transmembrane protein 54

Calculated MW:

24 kDa

Observed MW:

24 kDa, 43 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

TMEM54, also named as BCLP and CAC1, is widely expressed in cancer cell lines.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

