

For Research Use Only

# TUBGCP5 Polyclonal antibody, PBS Only

Catalog Number:14620-1-PBS



## Basic Information

Catalog Number:

14620-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6248

GenBank Accession Number:

BC071560

GeneID (NCBI):

114791

UNIPROT ID:

Q96RT8

Full Name:

tubulin, gamma complex associated protein 5

Calculated MW:

118 kDa

Observed MW:

118-120 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14620-1-AP (TUBGCP5 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14620-1-PBS in a different storage buffer formulation.



mouse brain tissue were subjected to SDS PAGE followed by western blot with 14620-1-AP (TUBGCP5 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14620-1-PBS in a different storage buffer formulation.