

For Research Use Only

# LPCAT3 Polyclonal antibody

Catalog Number: 14668-1-AP

3 Publications



## Basic Information

Catalog Number:

14668-1-AP

Concentration:

180 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6329

GenBank Accession Number:

BC065194

GeneID (NCBI):

10162

UNIPROT ID:

Q6P1A2

Full Name:

lysophosphatidylcholine  
acyltransferase 3

Calculated MW:

56 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Lijun Zou	39827610	Ecotoxicol Environ Saf	WB
Furui Feng	39260333	Biochem Biophys Res Commun	WB
Ling Kan Chi	38882385	Heliyon	WB, IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

